

■ MAB - INSERT

cod. MAB-ME-002



it is possible, using this insert, to adapt most common turbine tips to the various ultrasound machines on the market. Some possible uses of the MAB insert:

- a – in a big lift, to create antrostomy in sinus lift and elevate the Schneider membrane
- b – in a mini lift, to break into the maxillary sinus floor;
- c – preparation of the implant tunnel close to the lower alveolar nerve (use a full round head diamond point);
- d – in the finishing of the cavity surfaces in Conservative Dentistry. In the III class cavities, the use of a ball head diamond point divided in half (so as to avoid contact with the adjoining tooth) is advised;
- e – finishing of abutment preparation in fixed prosthesis (same point used in preparation, but with more fine-grained);
- f – resection of radicular apex in apicectomy;
- g – in endodontics, to remove the intra-cameral calcifications and to find the entrances of the radicular canals;
- h – to realize the microgeography and macrogeography on the vestibular surfaces for the fourth classes in the composite restorations of the frontal sector;
- i – separation of thin osseous crests (horizontal expansion of the crest).



■ MAB SPANNER

cod. MAB-03-001

